MicroDrop® **CalimeroJet**

MicroDrop® **MasterVent**

by **medel**

MPV TRUMA

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni per l'uso

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
ORDER DETAILS

Description

MicroDrop®CalimeroJet compressor-nebulizer system for children

MicroDrop®MasterVent compressor-nebulizer system for youth and adults

REF

M 51503-00

M 51403-00

Accessories

Compressor MicroDrop®CalimeroJet silver/red

M 50010-30000

Compressor MicroDrop®MasterVent silver/blue

M 50010-30100

CalimeroJet Nebulizer-Set consists of: CalimeroJet Nebulizer,

M 53700-00

Mouthpiece with valve, Tubing, child mask with adapter,

Air duct without valve, Stand, Instruction manual German/English

MasterVent Nebulizer-Set consists of: MasterVent Nebulizer,

M 53710-00

Mouthpiece with valve, Tubing, adult mask with adapter,

Air duct without valve, Stand, Instruction manual German/English

CalimeroJet Nebulizer without tubing

M 50020-33100

MasterVent Nebulizer without tubing

M 50020-33500

Year Set CalimeroJet Nebulizer consists of: 1xNebulizer CalimeroJet,

M 53700-03

1xTubing, 1xPorex Filter for MicroDrop compressor

Year Set MasterVent Nebulizer consists of: 1xNebulizer MasterVent,

M 53710-03

1xTubing, 1xPorex Filter for MicroDrop compressor

Mouthpiece CalimeroJet

M 50020-30700

Mouthpiece MasterVent

M 50020-30800

Medication cup CalimeroJet complete with nozzle

M 50020-33200

Medication cup MasterVent complete with nozzle

M 50020-33600

Nebulizer chamber yellow CalimeroJet

M 50020-30200

Nebulizer chamber blue MasterVent

M 50020-31400

Stand CalimeroJet

M 50020-30000

Stand MasterVent

M 50020-31200

Mask adult

M 50020-32000

Mask child

M 50020-32100

Air duct transparent without valve, for operation with mask

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Read the instructions

Alternating Current

CE 0123According to
93/42/EEC Directive. **IMPORTANT:**

Read these instructions carefully before using this device.

GENERAL INFORMATION

MicroDrop CalimeroJet/MasterVent is a compressor nebulizer system for aerosol therapy for intensive use. MicroDrop CalimeroJet/MasterVent is supplied with the CalimeroJet/MasterVent nebulizer, a new and patented exhalation/inhalation valves system, which increases the amount of inhaled medication, minimizing the medication waste.

INNOVATIVE DOUBLE-VALVE SYSTEM OPERATING DETAIL

During inhalation phase, the inhalation valve opens, automatically-operated by patient's breath, thereby allowing an additional air flow to the compressor-generated one, therefore increasing the nebulization of medication and the medication quantity effectively inhaled.

During exhalation phase, the inhalation valve closes avoiding medication to be dispersed in the environment while the exhalation valve opens thereby allowing the exhalation.

PRECAUTIONS

1. **Only use the device for its intended purpose as described in this instruction manual, therefore as a compressor nebulizer system for aerosoltherapy, following the indications of your doctor.** Any other form of use constitutes an improper use and is therefore dangerous. The manufacturer cannot be held liable for any damage caused by improper, incorrect or unreasonable use, or if the unit is connected to electrical installations which do not conform with safety regulations.
2. Keep this manual for future reference.
3. Do not operate the unit in presence of any anesthetic, inflammable mixtures or oxygen.
4. The correct device functioning can be affected by electromagnetic interference which exceeds the limits indicated by the European standards in force.
5. In case of failure and/or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" section. Do not handle nor open the compressor housing.
6. Any repairs must be carried out by an authorized dealer using only original spare parts. Safety of the device can be compromised in case the above mentioned indications are not met.
7. When using any electrical appliance certain important safety measures must always be observed, including the following:
 - use only original accessories and components;
 - never submerge the unit in water;
 - the device is not protected against water penetration;
 - never touch the unit when your hands are wet or moist;
 - do not leave the unit outdoor;
 - the unit, when operating, must be placed on a stable and horizontal surface;
 - make sure the air vent openings are not obstructed (pic. A/7);
 - do not allow unsupervised children or infirm persons to use the unit;
 - do not disconnect the unit by simply pulling out the plug from the wall socket.
8. Make sure that the electrical rating shown on the rating plate on the bottom of the unit corresponds to your main voltage and frequency before plugging in the device.
9. If the power plug provided with the device does not fit your wall socket, have the plug changed by an electrical appliance dealer. Do not use any adapter, simple or multiple, and/or extension cable. In case their use is necessary, make sure they are in compliance with safety standards, paying attention they do not exceed the maximum load.
10. Do not leave the unit plugged when not in use: unplug the device from the wall socket when not operating.
11. Follow the manufacturer's instructions for installing the device. The manufacturer cannot be held responsible for any damage caused to persons, animals or things by incorrect installation.
12. The power supply cord cannot be replaced by the user. In case the power supply cord becomes damaged, address to the technical servicing authorized by the manufacturer for replacement.
13. The power supply cord should always be fully unwound in order to prevent dangerous overheating.

14. Before performing any maintenance or cleaning operation, disconnect the device from the main supply, by unplugging it or by turning the main supply switch off.
15. Some parts of the unit are so small that they could be swallowed by children; do not let unattended children have access to the unit.
16. If you decide not to use the device any longer, it is recommended to dispose of it according to the regulations in force.
17. Remember to:
 - use this device only with medications prescribed by your doctor;
 - make the treatment only using the accessory recommended by your doctor for the therapy;

DEVICE SAFETY SYSTEMS

THERMAL PROTECTOR

The inhalation device is provided with a temperature protection device which is released when the device is overheating switching the unit off.

INSTRUCTIONS FOR USE

 **Remove the compressor from carrying case before use.**

Check the unit before each use, in order to check for possible function anomalies and/or damages due to transport and/or storage.

When inhaling, sit upright and relaxed at a table and not in armchair as otherwise you will constrict your respiratory passages and impair the effect of inhalation.

1. As soon as you have unpacked the device, check it for visible damages or defects. In particular, check for cracks in the plastic housing, which could expose electrical components.
2. Check for accessories integrity.
3. **Before using the device**, proceed with the cleaning operations as described in the "CLEANING AND DISINFECTION" section.
4. Open the nebulizer by turning the top in an anticlockwise direction (pic. D).
5. Make sure the water conduction cone is properly fitted onto the air conduction cone inside the nebulizer (pic. E).
6. Put the prescribed medication dose into the nebulizer (pic. F). The graduated scale can be useful to control the correct medication quantity.
7. Close the nebulizer by turning the two parts in a clockwise direction, paying attention that they are thoroughly sealed.
8. Connect one end of the air tube to the nebulizer (pic. G) and the other end to the air outlet on the device (pic. A/4).

9. Use the CalimeroJet/MasterVent nebulizer with inhalation valve only with mouthpiece having exhalation valve. When using the masks, it is not necessary to keep inhalation valve: remove it and connect the masks directly onto the nebulizer (pic. I).
10. Plug the device to the wall socket, making sure that the electrical rating shown on the rating plate on the bottom of the unit corresponds to your main voltage.
11. To start the treatment, switch the device on by pressing the switch (pic. A/2) in the front of the device onto the "I" position.
12. Inhale the aerosol solution using the prescribed accessory.
13. When the treatment has been completed, switch the unit off by setting the switch position "O" and disconnect the plug from the wall socket.
14. Wash the nebulizer and accessories as described in the "CLEANING AND DISINFECTION OF THE ACCESSORIES" section.
15. Put the nebulizer, accessories and power cord inside the carrying case (pic. A/14).

CLEANING AND MAINTENANCE

CLEANING OF THE MAIN UNIT

Use a clean dry cloth and mild detergents for cleaning the unit.

⚠ ATTENTION: during cleaning operations, make sure liquids do not get inside the device.

CLEANING AND DISINFECTION OF THE ACCESSORIES

After each treatment:

The nebulizer CalimeroJet/MasterVent, the mouthpiece must be sterilized by boiling in water (for max. 10 minutes) or in autoclave (max 121°C, see autoclave instructions for duration). We suggest you disinfect the masks using cold disinfecting liquids (solutions with hypochlorite). Dry the nebulizer using a clean cloth; after connect it to the air outlet (pic. A/4), switch the device on and let it work freely for 5-10 minutes.

⚠ Do not boil nor autoclave the air tube and the masks.

⚠ Following carefully the cleaning and disinfecting instructions is very important for the product performance and for a successful therapy.

REPLACEMENT OF THE NEBULIZER

The nebulizer must be replaced after a long period of inactivity, or when it is obstructed by dry medication, dust, etc.

Your nebulizer should last from 6 months to 1 year with general use.

Use original nebulizers only.

REPLACEMENT OF THE FILTER

The MicroDrop CalimeroJet/MasterVent Porex filter must be replaced approximately after 500 hours of use or each year; control the filter periodically (10-12 aerosol treatments); when it turns grey or brown or wet, replace it; open the filter holder (pic. A/5) and replace the filter with a new one. Don't try to clean the filter and to reuse.

Use original filters only.

⚠ Do not use the unit without filter.

POSSIBLE PROBLEMS AND HOW TO SOLVE THEM

If the device does not switch on:

1. make sure the plug is firmly fitted to the wall socket;
2. the thermal protector might have been released. The thermal protector may be released for the following reasons:
 - the device has been in operation for a too long time;
 - the device was used near heaters or in hot environments;
 - the air vent openings (pic. A/7) are obstructed and do not allow for the motor ventilation.

In this case, set the ON/OFF switch to "O" and let it rest for at least 15 minutes.

If the device does not nebulize, or only very weakly, check the following points:

1. make sure that the ends of the air tube are fitted tightly onto the main unit and onto the nebulizer;
2. check whether the nebulizer is empty from medication or full of an excessive medication quantity (MAX 7 ml);
3. check whether the nebulizer nozzle is obstructed;
4. check whether the air conduction cone is inserted (pic. E).

MAINTENANCE AND REPAIRS

In case of failure, address to qualified personnel authorized by MPV-TRUMA. **Never open the device.** There are no user-serviceable parts within. The compressor requires no lubrication or maintenance.

TECHNICAL SPECIFICATIONS

Model:	MicroDrop CalimeroJet/MasterVent
Electricity supply:	230 V, 50 Hz, 1A
Nebulization rate:	0,30 ml/min.
Particle size:	MMAD 1,9 µm (NaCl 0,9%) – API AEROSIZER MMD 4,90 µm – INSITEC/MALVERN
Nebulizer air flow:	3,3 l/min
Compressor air flow:	7,50 l/min
Compressor max. pressure:	250 kPa (2,50 bar)
Weight:	1,91 kg
Size:	(H) 106 x (L) 198 x (P) 233 mm
Noise Level:	max. 52 dBA

Class II Device (Double Safety Insulation)

Type B Device (device with specific protection against electrical hazards)

Device for continuous use.

Device not protected against sprinkles

Device not suitable for use in presence of anaesthetic substances inflammable with air, oxygen or nitrogen protoxide.

Device not suitable for use in anaesthesia or pulmonary ventilation systems.

Reference to Standards:	Electric Safety Standards CEI EN 60601.1 Electromagnetic Compatibility according to CEI EN 60601.1.2 Class II Medical Device according to 93/42/EEC "Medical Device" Directive.
Manufacturer:	Medel S.p.A., Via Micheli 7, 43056 S. Polo di Torrile (PR) Italy

Technical modifications reserved!

AMBIENT CONDITIONS

Storing Conditions

Temperature:	MIN -25°C - MAX 70°C
Humidity:	MIN 10%RH - MAX 95%RH

Operating Conditions

Temperature:	MIN 10°C - MAX 40°C
Humidity:	MIN 10%RH - MAX 95%RH
Altitude:	from 0 to 2000m above sea level

ORDER DETAILS

Description

MicroDrop®CalimeroJet compressor-nebulizer system for children	REF M 51503-00
MicroDrop®MasterVent compressor-nebulizer system for youth and adults	M 51403-00

Accessories

Compressor MicroDrop®CalimeroJet silver/red	M 50010-30000
Compressor MicroDrop®MasterVent silver/blue	M 50010-30100
CalimeroJet Nebulizer-Set consists of: CalimeroJet Nebulizer, Mouthpiece with valve, Tubing, child mask with adapter, Airduct without valve, Stand, Instruction manual German/English	M 53700-00
MasterVent Nebulizer-Set consists of: MasterVent Nebulizer, Mouthpiece with valve, Tubing, adult mask with adapter, Airduct without valve, Stand, Instruction manual German/English	M 53710-00
CalimeroJet Nebulizer without tubing	M 50020-33100
MasterVent Nebulizer without tubing	M 50020-33500
Year Set CalimeroJet Nebulizer consists of: 1xNebulizer CalimeroJet, 1xTubing, 1xPorex Filter for MicroDrop compressor	M 53700-03
Year Set MasterVent Nebulizer consists of: 1xNebulizer MasterVent, 1xTubing, 1xPorex Filter for MicroDrop compressor	M 53710-03
Mouthpiece CalimeroJet	M 50020-30700
Mouthpiece MasterVent	M 50020-30800
Medication cup CalimeroJet complete with nozzle	M 50020-33200
Medication cup MasterVent complete with nozzle	M 50020-33600
Nebulizer chamber yellow CalimeroJet	M 50020-30200
Nebulizer chamber blue MasterVent	M 50020-31400
Stand CalimeroJet	M 50020-30000
Stand MasterVent	M 50020-31200
Mask adult	M 50020-32000
Mask child	M 50020-32100
Airduct transparent without valvel, for operation with mask	M 50020-31900
Filter for MicroDrop compressor	M 50010-30200
Head cover CalimeroJet	M 50020-30400
Head cover MasterVent	M 50020-31500
Airduct CalimeroJet with valve	M 50020-33300
Airduct MasterVent with valve	M 50020-33700
Tubing CalimeroJet	M 50020-34000
Tubing MasterVent	M 50020-34100
Instruction manual MicroDrop CalimeroJet/MasterVent	M 50020-32900
Instruction manual CalimeroJet-/MasterVent-Vernebler	M 50020-33000

UK WARRANTY CERTIFICATE**Warranty conditions**

- The device is covered by a 5-year warranty from date of purchase for any defects in the design or in the materials used.
- The warranty consists in the replacement and/or repair free of charge of originally defective components.
- The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear.
- The device must be repaired by authorized technical service centres only.
- The device must be sent to the authorized service centre for repair within 8 days after the defect is noticed.
- The transport costs for the device are at user's charge.
- Any repair out of the warranty conditions are at user's charge.
- The warranty does not cover failures resulting from misuse, abuse of the product or in case the damage is not due to manufacturer's liability (accidental fall, rough transportation, etc.).
- The warranty does not involve any compensation for any kind of damages to persons or things, either direct or indirect, occurred when the device is not operating.
- The warranty is valid from the date of purchase certified by the ticket or the invoice which have to be peremptorily attached to the warranty coupon.
- The lack of warranty coupon properly filled in and validated by a certification of purchase voids the warranty.

**COUPON TO BE RETURNED IN CASE OF REPAIRS****THE WARRANTY IS VALID ONLY IF ACCOMPANIED BY THE TICKET/INVOICE**

Product type: _____

Model: _____

Serial No.: _____

Date of purchase: _____

BUYER DATA:

Name: _____

Address: _____

Phone number: _____

Defect description: _____

Signature for acceptance of the above mentioned warranty conditions

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